

## The Model DC2480 Knee CPM™



### *Extended Range of Motion*

The Model DC2480™ Knee CPM delivers an extended range of motion within a compact design that is ideal for any protocol. Boasting a range of -10 to 120 degrees and a sleek open carriage, the Model DC2480™ caters to the versatile needs of patients and providers alike.

Durable, mild steel construction allows the DC2480™ to accommodate home and hospital environments, while the easily understood control pendant provides reliable functionality. All Furniss Corp knee CPM models feature a unique space saving tripod position to reduce contamination and provide a convenient storage option. At just 27 inches long, the Model DC2480™ is proof that great things do come in small packages.

### Technical Specifications

Knee Flexion: 120 degrees  
Knee Extension: -10 degrees

Speed: 1 to 3 degrees/second  
Pause: 0 to 30 seconds

Limb Length: 21 to 41 inches  
Calf Length: 9 to 23 inches  
Thigh Length: 12 to 18 inches


Power Supply: 110v  
Unit Weight: 25 lbs.  
Ship Weight: 38 lbs.  
Length: 27 inches

### **CPM Devices + Softgoods**

*Manufactured for Providers All Over the World*

US 800.273.5233 Fax 614.871.9353

KneeCPM.com



# DC2480 Knee CPM

## Pediatric Athletic Bariatric

Proper anatomical therapy, accurate range of motion and ultimate comfort for all patient types

- Lightweight (25 lbs.) and Portable (28 in.)
- Constructed with mild steel and a dual ball screw assembly—exclusively
- Circuitry located in base of unit instead of hand control for fundamental protection
- Easy to understand hand-held control features a bold LED digital display
- Patient lock out feature prevents inadvertent changes to settings
- Compliance meter records treatment time
- Features safety stop if motion is obstructed
- Universal left/right design
- Quick conversion from adult to pediatric—no additional accessories needed
- Tripod position offers easy storage option and reduces contamination
- Axis to axis motion, skid pad and bed stabilizing system prevent migration
- Knee CPM softgoods are constructed with deep pile synthetic sheepskin and high tensile strength Velcro to ensure patient comfort and proper positioning during therapy. (order #1828-1)
- Two-year extended manufacturer warranty
- HCPCS E0935



### Indications

Total knee replacement  
ACL reconstruction  
Manipulation under anesthesia  
Ligament reconstruction  
Meniscus repair  
Prosthetic replacement  
Prolonged joint immobilization

Class I Medical Device

### Furniss Corporation Defined

As a US based medical manufacturer, The Furniss Corporation builds Continuous Passive Motion (CPM) devices for the upper and lower extremities. The durable, user friendly devices are built to withstand the rigorous hospital and home care environment. Furniss Corp also manufactures patient softgoods and offers a comprehensive technical service and repair center. Expect same-day shipping, immediate access to educated product professionals and strong product lines with exceptional performance. Visit the website at [KneeCPM.com](http://KneeCPM.com) to experience the interactive tour.

CPM Devices + Softgoods  
US 800.273.5233  
[KneeCPM.com](http://KneeCPM.com)